Product name - VACUUM MATTRESS PUMP
Cod: POM002 / POM003
Regulations: The products defined as medical devices are compliant with MDR 745/2017 of 05/04/2017





Italiano

Descrizione
Pompa manuale per materasso a depressione, disponibile in versione con corpo in PVC o in alluminio verniciato. Completa di tubo con raccordo e manico girevole, utilizzabile in qualsiasi posizione. La versione in metallo è completamente smontabile.

Dimensioni
POM002: dimensioni 110 x 220 x 450 mm; peso 1,05 kg; portata 2 x 2000 cc; pressione massima 500 mbar (7.2 psi); temperatura di funzionamento -10°C / +40°C.
POM003: lunghezza totale 485 mm; lunghezza camera 290 mm; diametro camera 60 mm; lunghezza tubo 800 mm; peso 1,5 kg; capacità 0,75 litri.

Istruzioni d'uso
Collegare il tubo al materasso a depressione mediante il raccordo. Azionare la pompa tramite il manico girevole fino al raggiungimento del vuoto desiderato. Utilizzare il dispositivo su superficie stabile e verificare il corretto collegamento del tubo prima dell'uso. Al termine, riporre in luogo asciutto e pulito. Per la versione POM003, se necessario, smontare i componenti per agevolare trasporto e manutenzione.
(Istruzioni d'uso formulate su base tecnica coerente con la destinazione del prodotto.)

English

Description
Manual pump for vacuum mattress, available with PVC body or painted aluminium body. Supplied with hose and fitting and rotating handle, operable in any position. The metal version is fully dismountable.

Dimensions
POM002: dimensions 110 x 220 x 450 mm; weight 1.05 kg; capacity 2 x 2000 cc; max pressure 500 mbar (7.2 psi); working temperature -10°C / +40°C.
POM003: total length 485 mm; cylinder length 290 mm; cylinder diameter 60 mm; hosepipe length 800 mm; weight 1.5 kg; capacity 0.75 litre.

Instructions for use
Connect the hose to the vacuum mattress through the fitting. Operate the pump using the rotating handle until the required vacuum level is reached. Use the device on a stable surface and check the correct hose connection before use. After use, store in a clean, dry place. For POM003, components may be disassembled if required for easier transport and maintenance.
(Instructions for use drafted in line with the intended use of the product.)

Français

Description
Pompe manuelle pour matelas à dépression, disponible avec corps en PVC ou en